

# EASYSHOE VERSA APPLICATION

**The EasyShoe Versa is a hybrid nail-on shoe that features the concussion absorption powers of polyurethane and an internal steel insert. The steel insert gives the shoe more stability in the front portion and offers a stopping point when driving nails.**

The steel insert is located in the toe area and continues to the widest part of the hoof. This positioning allows the rear portion of the hoof the ability to flex with each footfall. The Versa offers incredible versatility for distance rides, trail horses, and horses working on hard surfaces. The wide web design also makes the Versa ideal for therapy cases because it holds pour-in packing very well.

The EasyShoe Versa can be narrowed with a hammer or opened with a shoe spreader. Although the shoe can be shaped, we anticipate most farriers will simply drop it on and shave away excess urethane material. There are only 4mm between each shoe size, which makes urethane removal quick and easy. The goal is to fit each hoof with little to no modification required to save farriers time and money.

No predrilling needed on 5 of the 6 nail locations on the Versa. 5 of the 6 nailing locations are open and for nailing. We did sneak an extra nail location in behind the clip that will need to be drilled if you choose to nail in that location.

*Sold in Pairs.*

## EasyShoe Versa Installation Instructions

The EasyShoe Versa has an internal steel core. The EasyShoe Versa can be narrowed with a hammer or opened with a shoe spreader. Although the shoe can be shaped, we anticipate most farriers will simply drop it on and shave away excess urethane material. There are only 4mm between each shoe size, which makes urethane removal quick and easy. The goal is to fit each hoof with little to no modification required to save farriers time and money.



1.



2.



5.

1. Place the shoe on the hoof and trace with a marker.
2. Remove excess material. A belt sander will work well. A band saw can also be used followed by a belt sander for a smooth finish.
3. Both metal and urethane can be removed in the fitting process. Bend side clips to match hoof angle.
4. Make adjustments to heel length if needed.
5. We recommend E4 Slim for Versa sizes (114-134) E5 Slim sizes for Versa sizes (138-154). Make sure to seat nail heads against steel in shoe. Finish and clinch in the usual fashion.



**EasyCare Inc.**  
The Ultimate in Hoof Protection

# HOW TO MEASURE FOR THE EASYSHOE VERSA

# EasyShoe **VERSA**

## APPLICATION GUIDE

### The EasyShoe Versa features:



A hind pattern that fits.



Front and hind nailing patterns are sized to fit with limited modifications.



Stronger, sleeker, and more traditional clips.



More nailing options with bonus nailing location behind the clip if needed.



High abrasion.



Easily adjust breakerover.

- Trail Riding
- Long Distance Events
- Speed Events
- Driving/Parades/Performance

### Join the Conversation

For tips and tricks follow EasyCare:




### A product of EasyCare Inc.

EasyCare Guarantee

800.447.8836 | easycareinc.com

All products should be used in accordance with the instructions in this application guide. EasyCare Inc. hereby disclaims responsibility for the misuse of any product and for losses incurred through the use of any product described in this guide, on EasyCare Inc. product packaging, or on easycareinc.com. EasyCare Inc. products are not intended to diagnose or cure any medical condition.

1. After a fresh trim, measure the width of the hoof across the sole at the widest point.
2. Measure the length of the hoof from the toe to the buttress line of the heel. The buttress line is the farthest weight bearing point of the heel where the hoof wall ends. Do not include the heel bulbs in the measurement.
3. Compare your measurement with the size chart below.



4. Ideally, the length and width measurement will fit into the same size. If the length and width measurements indicate different sizes, select the larger size.

Size	Width	Length	Thickness	Stock Status
102	96-102mm			Coming Soon
106	100-106mm			Coming Soon
110	104-110mm			Coming Soon
114	108-114mm			Coming Soon
118	112-118mm	127mm	14.5mm	In Stock
122	116-122mm	132mm	14.5mm	In Stock
126	120-126mm	136mm	14.5mm	In Stock
130	124-130mm	140mm	14.5mm	In Stock
134	128-134mm	145mm	15mm	In Stock
138	132-138mm	149mm	15mm	In Stock
142	136-142mm	153mm	16mm	In Stock
146	140-146mm	158mm	16mm	In Stock
150	144-150mm	162mm	16mm	In Stock
154	148-154mm	166mm	16mm	In Stock
158	152-158mm			Coming Soon
102	96-102mm			Coming Soon
106	100-106mm			Coming Soon
110	104-110mm			Coming Soon
114	108-114mm	130mm	14.5mm	In Stock
118	112-118mm	134mm	14.5mm	In Stock
122	116-122mm	139mm	14.5mm	In Stock
126	120-126mm	143mm	14.5mm	In Stock
130	124-130mm	148mm	14.5mm	In Stock
134	128-134mm	152mm	14.5mm	In Stock
138	132-138mm	157mm	15mm	In Stock
142	136-142mm	161mm	15mm	In Stock
146	140-146mm	166mm	16mm	In Stock
150	144-150mm	171mm	16mm	In Stock
154	148-154mm			Coming Soon
158	152-158mm			Coming Soon



**EasyCare Inc.**  
The Ultimate in Hoof Protection