


Easyboot Edge Application



****Note:** The Easyboot Edge is designed for a barefoot hoof that is slightly longer than wide and has limited flare. The Easyboot Edge does not work well on a round hoof or a hoof with excess flare.

- 1. Loosen the Easyboot Edge worm clamp** using the EasyKey by turning it counter-clockwise (see illustration).
- 2. Fold the gaiter back before fitting the boot over the hoof.**
- 3. Pick up the horse's leg and support the cannon bone between your knees** so that you can use two hands when applying the boot.
- 4. Using technique rather than strength, place the boot over the horse's toe** and, using a slight back and forth twisting motion, slide the boot on the rest of the way. You can also put the foot down on the ground and allow the horse to put weight into it. Don't use the gaiter to pull the boot on. If you are having trouble please refer to the fitting video found at easycareinc.com.




- 5. Check for proper fit before tightening worm clamp.** Make sure all portions of the boot are below the hairline. If all portions of the boot (excluding the gaiter) are below the hairline, tighten the worm clamp with the EasyKey by turning it clockwise (see illustration). Snug the boot, but **don't overtighten**. The worm clamp system is designed so that a flathead screwdriver or even a penny can be used to tighten or loosen the worm clamp if the EasyKey is misplaced.



- 6. Fasten the Easyboot Edge Gaiter.** Pull the top of the gaiter tight around the pastern. Fasten the hook and loop straps. The tighter the fit, the better the boot will perform.

Due to their great deal of elasticity, gaiter straps cannot be over tightened. Overlap the hook and loop straps as much as needed to insure a tight fit and avoid rubbing issues. This overlapping does not compromise the integrity of the gaiter. When properly tightened, it should be difficult to fit a finger between the pastern and the top of the gaiter.



Removing the Easyboot Edge.

- Loosen the worm clamp using the EasyKey by turning it counter-clockwise (see illustration).
- Unfasten the gaiter straps.
- Slip the boot off the hoof.



Have Questions? Need Help? Visit easycareinc.com



EasyCare
easycareinc.com

EasyCare Inc. Warranty

EasyCare offers a 90 day repair or replacement warranty for all hoof boots (purchased as new) manufactured under the EasyCare label. Please register your product at http://www.easycareinc.com/warranty_reg.aspx and keep your original dated receipt(s) as proof of purchase during the warranty period.

If a boot fails during warranty, it is generally a minor hardware issue. Please bear in mind that most components of EasyCare hoof boots are either replaceable or repairable. This includes wire laces, cables, buckles, gears, rollers, screws and straps. Most repairs can be easily accomplished with simple household tools. Depending on the circumstances, EasyCare will send you repair parts, or ask that you return the boots to us for evaluation. Videos of many repairs are available on our website at www.easycareinc.com.

While this warranty is very broad in nature, it does not cover wear (including sole breakthrough when worn over horseshoes), misuse, abuse, factory seconds, bargain bin items, or boots purchased as used. Products returned for evaluation must be clean, dry and free of debris, dirt, sand and manure. *Failure to return clean products will result in product being shipped back at sender's expense or a cleaning fee may be assessed.* Damage inflicted by omnivorous horses, rodents or other unusual abuse is not covered.

Returns must be authorized before sending and must be returned with an RMA number. Contact EasyCare directly for all warranty issues at 1-800-447-8836 to receive an RMA number.

Specifically:

1. Returns for credit only (refund of purchase price only): Products must be in new, unused condition, free of dirt and wear marks. It is your responsibility to evaluate the fit of your new product in a manner where you may return it in new condition if you have doubts about the fit. Be certain that you know the return policy of product purchased from dealers other than direct from EasyCare (the manufacturer). EasyCare will not give credit for product purchased elsewhere. Return to point of purchase in that circumstance. Also, all brochures, inserts

Online Warranty Registration:
easycareinc.com/warranty_reg.aspx

and/or free gifts must accompany returned items.


2. Exchanges of new product (size or color issue): EasyCare will only exchange NEW items with same product style and model for same product when requested with NEW items only. You may exchange for a different size or color exchange only, no exceptions. If you wish to have a different style of hoof boot, etc., you must return new item for refund (original point of purchase) and then purchase new choice.
3. Warranty failure (used product that is within 90 days of purchase): Contact us at 800-447-8836. EasyCare will either provide repair parts directly to you, or repair/replace after our evaluation. Items must be returned to EasyCare at the consumers cost with original receipt. EasyCare will return the items via ground service (UPS or FEDX) free of charge. Requests for faster shipping service will be at the expense of the consumer.
4. Not under warranty (over 90 days or no proof of purchase): Repair parts may be purchased from your local dealer or EasyCare, or product may be returned for evaluation and possible repair (parts, labor and freight charges will apply).
5. Call for an RMA number (under warranty or not): Returns must be authorized before sending and must be returned with an RMA number clearly written on outside of packaging. Contact us at 800-447-8836 to receive an RMA number. Returns will be refused without an RMA number on outside of package.

6. Freight charges (outside warranty period): Freight charges to and from EasyCare for repairs, exchanges and items outside the warranty period will be the purchaser's responsibility.
7. EasyCare staff is always willing to discuss any specific issues you may have regarding the performance of your EasyCare or Old Mac's products. Please contact us at 800-447-8836 or email to admin@easycareinc.com. Visit our website, easycareinc.com, for tips and product information that will help with the use of any of our product line.




How to Measure for Boots

1. After a fresh trim, measure the width of the hoof across the bottom at the widest point.
2. Measure the length of the hoof from the toe to the buttress line of the heel. The buttress line is the farthest weight bearing point of the heel where the hoof wall ends. Do not include the heel bulbs in the measurement.



Measuring for
Hoof Width



Measuring for
Hoof Length

3. Compare your measurement with the size chart below. Please note that while most of our other boots have their own unique size charts, all of the new 2009 boots all use the same size chart.
4. Ideally, the length and width measurement will fit into the same size. If the length and width measurements indicate different sizes, select the larger size. If the width and length measurements are different by more than one size, your selected boot style is not recommended. Be sure to check your measurements against an alternative boot in the EasyCare line—another style may accommodate your horse's hoof shape better.

| Size | Width in inches | Length in inches | Size | Width in mm | Length in mm |
|------|------------------|--------------------|------|-------------|--------------|
| 0 | 4 3/8" - 4 9/16" | 4 5/8" - 4 15/16" | 0 | 111-116mm | 118-124mm |
| 1 | 4 9/16" - 4 7/8" | 4 15/16" - 5 1/8" | 1 | 117-123mm | 125-130mm |
| 2 | 4 7/8" - 5 1/8" | 5 1/8" - 5 9/16" | 2 | 124-130mm | 131-140mm |
| 3 | 5 1/8" - 5 1/2" | 5 9/16" - 5 15/16" | 3 | 131-139mm | 141-150mm |

We continue to strive for the best products on the market with you and your horse in mind. We hope you enjoy our product selection and website improvements. Thank you for your continued support and patronage.

The Fine Print:

EasyCare Inc. hereby disclaims responsibility for the misuse of any products described in this brochure and web site and/or losses incurred through the use of products by any method other than that which is outlined in this catalog or web site or EasyCare's product packaging. All products should be used in accordance with the instructions printed on the product's label. EasyCare Inc products are not intended to treat, cure or diagnose any medical condition.

easycareinc.com 800.447.8836



Easyboot® Edge

Application Guide