


Easyboot Glove and Glue-On Hoof Boot Measurement Log

Keep track of your horse's hoof growth.
 Questions? Give us a call! 1.800.447.8836.



Measure the width of the hoof across the bottom at the widest point.

Do not assume that, if you wear a certain size in another boot style, you will wear the same size Glue-On or Glove.

Success of the Easyboot Glove and Easyboot Glue-On are determined by accurate measurements and a snug fit.

We have found that measurements are more accurate and easier to take using the metric system.

Size	Width in mm	Length in mm
00	96-99mm	106-109mm
00.5	100-103mm	110-113mm
0	104-107mm	114-117mm
0.5	108-111mm	118-121mm
1	112-115mm	122-125mm
1.5	116-119mm	126-129mm
2	120-123mm	130-133mm
2.5	124-127mm	134-137mm
3	128-131mm	138-141mm
3.5	132-135mm	142-145mm
4	136-139mm	146-149mm



Measure the length of the hoof from the toe to the buttress line of the heel. The buttress line is the farthest weight bearing point of the heel where the hoof wall ends. Do not include the heel bulbs in the measurement.

HORSE'S NAME:

Date: _____

	Width	Length	Boot Size
Right Front			
Left Front			
Right Hind			
Left Hind			

Date: _____

	Width	Length	Boot Size
Right Front			
Left Front			
Right Hind			
Left Hind			

Date: _____

	Width	Length	Boot Size
Right Front			
Left Front			
Right Hind			
Left Hind			

Date: _____

	Width	Length	Boot Size
Right Front			
Left Front			
Right Hind			
Left Hind			



EasyCare